

## SCIALYTIQUE PLAFONNIER LED700-500



**Star**Care  
M E D I C A L

### TECHNICAL SPECIFICATIONS FOR LED LIGHT

| Technical specifications for lights        | LED700             | LED500             |
|--------------------------------------------|--------------------|--------------------|
| Central illuminance at 1m working distance | 40,000~16,0000 lux | 40,000~130,000 lux |
| Focusable light-field size at 1 m distance | 120~300mm          | 120~300mm          |
| Colour temperature                         | 3500~5000K         | 3500~5000K         |
| Depth of Focus                             | >700mm             | >700mm             |
| Colour Index Ra                            | 95                 | 95                 |
| Power Supply                               | 110~240V 50/60Hz   | 110~240V 50/60Hz   |
| Effective bulb service life                | >50,000h           | >50,000h           |
| Diameter of lamp housing                   | 700mm              | 500mm              |
| Power consumption of light heads(W)        | 60                 | 48                 |
| Installation height (Ceiling type)         | 280-320 cm         | 280-320 cm         |